

Anti-CD35 antibody [E11] (STJ16101912)

STJ16101912

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD35 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/FC/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	E11
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human acute monocytic leukaemia cells
Immunogen Region	
Specificity	
Immunogen Sequence	