

Anti-CD45RA antibody [F8-11-13] (STJ16101916)

STJ16101916

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD45RA is suitable for use in Immunohistochemistry, Flow Cytometry and Immunohistochemistry research applications. |
| Applications | IHC-F/FC/IHC-P |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | F8-11-13 |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified |
| Dilution Range | |
| Formulation | PBS with <0.1% sodium azide |
| Isotype | IgG1 |
| Storage | Store aliquots at -20°C, and thawed aliquots at 2-8°C. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|---------------------|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Human T lymphocytes |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |