

Anti-CD72 antibody [Bu40] (STJ16101925)

STJ16101925

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD72 is suitable for use in Immunohistochemistry, Flow Cytometry, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/FC/IP/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Bu40
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	971
Gene Symbol	CD72
Uniprot ID	CD72_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	