

Anti-Chymotrypsin antibody [CHYMO10-104.5] (STJ16101930)

STJ16101930

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Chymotrypsin is suitable for use in Immunohistochemistry, ELISA, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/ELISA/IF/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	CHYMO10-104.5
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG3
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified human pancreatic chymotrypsin
Immunogen Region	
Specificity	
Immunogen Sequence	