

Anti-HIF1 Alpha antibody [Halpha111a] (STJ16101933)

STJ16101933

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HIF1 Alpha is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Halpha111a
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant human HIF1 alpha
Immunogen Region	
Specificity	
Immunogen Sequence	