

Anti-Lysozyme antibody [BGN/0696/5B1] (STJ16101935)

STJ16101935

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lysozyme is suitable for use in ELISA, Flow Cytometry and Immunohistochemistry research applications.
Applications	ELISA/FC/IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BGN/0696/5B1
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Native, from urine
Immunogen Region	
Specificity	
Immunogen	
Sequence	