

Anti-MAdCAM-1 antibody [17F5] (STJ16101936)

STJ16101936

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MAdCAM-1 is suitable for use in ELISA, Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	ELISA/FC/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Camel

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17F5
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant soluble MAdCAM-1-fc fusion protein
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	