

Anti-N-Cadherin antibody [13A9] (STJ16101940)

STJ16101940

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-N-Cadherin is suitable for use in Immunohistochemistry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/IF/IP/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	13A9
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant MBP fusion protein containing the entire cytoplasmic domain of human N-cadherin
Immunogen Region	
Specificity	
Immunogen Sequence	