

Anti-Parathyroid Hormone antibody [BGN/1F8] (STJ16101942)

STJ16101942

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Parathyroid Hormone is suitable for use in Immunohistochemistry, ELISA and Immunohistochemistry research applications.
Applications	IHC-F/ELISA/IHC-P
Host/Source	Mouse
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BGN/1F8
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids 1-34 of mature PTH conjugated to a proprietary carrier molecule
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	