

Anti-Protein Gene Product 9.5 antibody [31A3] (STJ16101944)
STJ16101944**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Protein Gene Product 9.5 is suitable for use in ELISA, Immunohistochemistry and Western Blot research applications.
Applications	ELISA/IHC-P/WB
Host/Source	Mouse
Reactivity	Human/Rat/Rabbit/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	31A3
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Native, from brain
Immunogen Region	
Specificity	
Immunogen Sequence	