

Anti-TALIN-1 antibody [97H6] (STJ16101946)

STJ16101946

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-TALIN-1 is suitable for use in Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IF/IP/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	97H6
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	his-tagged polypeptide containing residues 489–911 of mouse talin1
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	