

**Anti-TALIN-2 antibody [53.8] (STJ16101947)**

STJ16101947

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-TALIN-2 is suitable for use in Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IF/IP/IHC-P/WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Rabbit

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	53.8
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Purified
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with <0.1% sodium azide
<b>Isotype</b>	IgG2b
<b>Storage</b>	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Recombinant Mouse talin-2 (aa2300-2542).
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	