

Anti-Tenascin antibody [T2H5] (STJ16101948)

STJ16101948

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Tenascin is suitable for use in Immunohistochemistry, ELISA, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F/ELISA/IP/IHC-P/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	T2H5
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Dilution Range	
Formulation	PBS with <0.1% sodium azide
Isotype	IgG1
Storage	Store aliquots at -20°C, and thawed aliquots at 2-8°C.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified protein from human mammary tumour
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	