

Anti-ASGPR/CD301 antibody [125A10.03] {Biotin} (STJ16101990)
STJ16101990**GENERAL INFORMATION**

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-ASGPR/CD301 is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC/WB |
| Host/Source | Mouse |
| Reactivity | Human/Dog |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 125A10.03 |
| Concentration | 0.5 mg/mL |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |