

Anti-DEC205/CD 205 (NLDC145) antibody [NLDC145] (STJ16102000)

STJ16102000

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-DEC205/CD 205 (NLDC145) is suitable for use in FACS, Immunohistochemistry and Immunofluorescence research applications.
Applications	FACS/IHC-F/IF
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	NLDC145
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence