

Anti-CD303 antibody [102G7] {Alexa-fluor-488} (STJ16102026)
STJ16102026**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD303 is suitable for use in Flow Cytometry and Immunofluorescence research applications.
Applications	FC/IF
Host/Source	Mouse
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	102G7
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence