

Anti-Birch Pollen antibody [110G10.01] (STJ16102037)

STJ16102037

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Birch Pollen is suitable for use in Western Blot and Immunoprecipitation research applications.
Applications	WB/IP
Host/Source	Mouse
Reactivity	Pollen

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	110G10.01
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	