

Anti-Beta 2-microglobuline antibody [mAb4] (STJ16102039)

STJ16102039

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Beta 2-microglobuline is suitable for use in FACS and ELISA research applications.
Applications	FACS/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	mAb4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	