

Anti-Bet V1 antibody [H.G28C10.01] (STJ16102043)
STJ16102043**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-Bet V1 is suitable for use in ELISA and RImmunoprecipitationA research applications.
Applications	ELISA/RIPA
Host/Source	Human
Reactivity	Pollen

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H.G28C10.01
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG4k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	