

Anti-Bet V1 antibody [H.G5F2A4.09] (STJ16102044)
STJ16102044**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-Bet V1 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Human
Reactivity	Pollen

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H.G5F2A4.09
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence