

Anti-CD1a antibody [214A9.01] {Alexa-fluor-546} (STJ16102048)

STJ16102048

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD1a is suitable for use in FACS research applications.
Applications	FACS
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	214A9.01
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgM
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	909
Gene Symbol	CD1A
Uniprot ID	CD1A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	