

Anti-CD1a antibody [201B5.08] {Alexa-fluor-488} (STJ16102052)
STJ16102052**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD1a is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Mouse
Reactivity	Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	201B5.08
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	909
Gene Symbol	CD1A
Uniprot ID	CD1A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	