

Anti-CD14 antibody [MOP9.25] (STJ16102058)

STJ16102058

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD14 is suitable for use in FACS, Cell Universal and Immunohistochemistry research applications.
Applications	FACS/Cell Depletion/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MOP9.25
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b, k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	929
Gene Symbol	CD14
Uniprot ID	CD14_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	