

Anti-CD54 antibody [117G12.02] {Alexa-fluor-488} (STJ16102067)

STJ16102067

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD54 is suitable for use in Flow Cytometry, Immunofluorescence and ELISA research applications.
Applications	FC/IF/ELISA
Host/Source	Mouse
Reactivity	Human/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	117G12.02
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence