

**Anti-CD54 antibody [117G12.02] {Alexa-fluor-647} (STJ16102069)**

STJ16102069

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD54 is suitable for use in Flow Cytometry, Immunofluorescence and ELISA research applications.
<b>Applications</b>	FC/IF/ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Sheep

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	117G12.02
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Alexa-fluor-647
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20°C for up to 1 year.

**TARGET INFORMATION**

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>