

**Anti-CD54 antibody [217E2.10] (STJ16102074)**

STJ16102074

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-CD54 is suitable for use in FACS, Immunofluorescence, Immunohistochemistry and Immunohistochemistry research applications.
<b>Applications</b>	FACS/IF/IHC-F/IHC-P
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	217E2.10
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at -20°C for up to 1 year.
<b>Instruction</b>	

**TARGET INFORMATION**

**Gene ID**  
**Gene Symbol**  
**Uniprot ID**  
**Immunogen**  
**Immunogen**  
**Region**  
**Specificity**  
**Immunogen**  
**Sequence**