

Anti-DCIR/CD367 antibody [108H8.03] (STJ16102085)
STJ16102085**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DCIR/CD367 is suitable for use in FACS and ELISA research applications.
Applications	FACS/ELISA
Host/Source	Mouse
Reactivity	Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	108H8.03
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	