

Anti-DC-LAMP/CD208 antibody [104G4] {Alexa-fluor-488} (STJ16102086)

STJ16102086

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DC-LAMP/CD208 is suitable for use in Immunofluorescence and FACS research applications.
Applications	IF/FACS
Host/Source	Mouse
Reactivity	Human/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	104G4
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence