

Anti-DC-LAMP/CD208 antibody [104G4] {Biotin} (STJ16102089)
STJ16102089**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DC-LAMP/CD208 is suitable for use in FACS and Immunohistochemistry research applications.
Applications	FACS/IHC-F
Host/Source	Mouse
Reactivity	Human/Sheep

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	104G4
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	