

Anti-DC-LAMP/CD208 antibody [1010E1.01] {Biotin} (STJ16102094)
STJ16102094**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-DC-LAMP/CD208 is suitable for use in FACS and Immunohistochemistry research applications.
Applications	FACS/IHC-F
Host/Source	Rat
Reactivity	Mouse/Dog/Cat/Human/Sheep/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1010E1.01
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence