

Anti-DC-LAMP/CD208 antibody [208B5] {Alexa-fluor-546} (STJ16102107)

STJ16102107

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-DC-LAMP/CD208 is suitable for use in FACS research applications. |
| Applications | FACS |
| Host/Source | Mouse |
| Reactivity | Mouse/Human |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 208B5 |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-546 |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG3 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |