

Anti-DC-SIGN/CD209 antibody [102E11.06] {Biotin} (STJ16102120)
STJ16102120**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DC-SIGN/CD209 is suitable for use in FACS, Immunohistochemistry and Immunohistochemistry research applications.
Applications	FACS/IHC-F/IHC-P
Host/Source	Mouse
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	102E11.06
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2b, k
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence