

Anti-DC-SIGN/CD209 antibody [109H12.03] (STJ16102140)

STJ16102140

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DC-SIGN/CD209 is suitable for use in Flow Cytometry, Immunohistochemistry and Immunofluorescence research applications.
Applications	FC/IHC/IF
Host/Source	Mouse
Reactivity	Human/Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	109H12.03
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence