

Anti-DERP-1 antibody [H.AK6.A3.F4] (STJ16102149)

STJ16102149

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-DERP-1 is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA/WB
Host/Source	Human
Reactivity	House Dust Mite

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H.AK6.A3.F4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence