

Anti-DERP-1 antibody [H.AA.IV7A3] {Biotin} (STJ16102150)
STJ16102150**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-DERP-1 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Human
Reactivity	House Dust Mite

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H.AA.IV7A3
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	