

**Anti-DERP-1 antibody [H.AA.IV7A3] {HRPO} (STJ16102151)**  
STJ16102151**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Human monoclonal antibody anti-DERP-1 is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Human
<b>Reactivity</b>	House Dust Mite

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	H.AA.IV7A3
<b>Concentration</b>	0.5 mg/mL
<b>Conjugation</b>	HRPO
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	
<b>Formulation</b>	TRIS, PBS, pH7 and antibiotics
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Store at -20°C for up to 1 year.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	