

Anti-DORA/IgSF6 antibody [104A10.01] {Alexa-fluor-647} (STJ16102161)

STJ16102161

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DORA/IgSF6 is suitable for use in FACS and Western Blot research applications.
Applications	FACS/WB
Host/Source	Mouse
Reactivity	Human/Mouse/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	104A10.01
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence