

**Anti-MMP-12/HME antibody [603E6.22] (STJ16102176)**

STJ16102176

**GENERAL INFORMATION**

|                          |                                                                                                                        |
|--------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                     |
| <b>Short Description</b> | Mouse monoclonal antibody anti-MMP-12/HME is suitable for use in Immunohistochemistry and ELISA research applications. |
| <b>Applications</b>      | IHC-F/ELISA                                                                                                            |
| <b>Host/Source</b>       | Mouse                                                                                                                  |
| <b>Reactivity</b>        | Human                                                                                                                  |

**PRODUCT PROPERTIES**

|                            |                                  |
|----------------------------|----------------------------------|
| <b>Clonality</b>           | Monoclonal                       |
| <b>Clone ID</b>            | 603E6.22                         |
| <b>Concentration</b>       | 0.5 mg/mL                        |
| <b>Conjugation</b>         | Unconjugated                     |
| <b>Purification</b>        | Affinity Purified                |
| <b>Dilution Range</b>      |                                  |
| <b>Formulation</b>         | TRIS, PBS, pH7 and antibiotics   |
| <b>Isotype</b>             | IgG1                             |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year. |

**TARGET INFORMATION**

|                         |
|-------------------------|
| <b>Gene ID</b>          |
| <b>Gene Symbol</b>      |
| <b>Uniprot ID</b>       |
| <b>Immunogen</b>        |
| <b>Immunogen Region</b> |
| <b>Specificity</b>      |
| <b>Immunogen</b>        |
| <b>Sequence</b>         |

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081