

Anti-IL3R Alpha/CD123 antibody [107D2.08] {Alexa-fluor-488} (STJ16102186)

STJ16102186

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IL3R Alpha/CD123 is suitable for use in FACS and Immunofluorescence research applications.
Applications	FACS/IF
Host/Source	Mouse
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	107D2.08
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	