

Anti-IL3R Alpha/CD123 antibody [107D2.08] {Alexa-fluor-546} (STJ16102187)

STJ16102187

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IL3R Alpha/CD123 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	107D2.08
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence