

Anti-IL3R Alpha/CD123 antibody [107D2.08] {Alexa-fluor-647} (STJ16102188)

STJ16102188

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-IL3R Alpha/CD123 is suitable for use in FACS and Immunofluorescence research applications. |
| Applications | FACS/IF |
| Host/Source | Mouse |
| Reactivity | Human/Dog |

PRODUCT PROPERTIES

| | |
|----------------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 107D2.08 |
| Concentration | 0.5 mg/mL |
| Conjugation | Alexa-fluor-647 |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage Instruction | Store at -20°C for up to 1 year. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |