

Anti-IL3R Alpha/CD123 antibody [107D2.08] (STJ16102190)

STJ16102190

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IL3R Alpha/CD123 is suitable for use in Immunohistochemistry, Western Blot, FACS and Immunoprecipitation research applications.
Applications	IHC-P/WB/FACS/IP
Host/Source	Mouse
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	107D2.08
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence