

Anti-CD39 antibody [AC2.5] (STJ16102205)

STJ16102205

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD39 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	FC/IHC-F/IP/WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AC2.5
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage	Store at -20°C for up to 1 year.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence