

Anti-IRAK-4 antibody [7D8] {Alexa-fluor-546} (STJ16102226)

STJ16102226

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IRAK-4 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	7D8
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence