

Anti-HKCR antibody [117E10.10] {Alexa-fluor-647} (STJ16102234)

STJ16102234

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HKCR is suitable for use in FACS and Immunohistochemistry research applications.
Applications	FACS/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	117E10.10
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	