

Anti-Langerin/CD207 antibody [808E10.01] {Alexa-fluor-488} (STJ16102236)
STJ16102236**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Langerin/CD207 is suitable for use in FACS research applications.
Applications	FACS
Host/Source	Mouse
Reactivity	Human/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	808E10.01
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-488
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	