

Anti-Langerin/CD207 antibody [310F7.02/HD26] {Alexa-fluor-546} (STJ16102242)

STJ16102242

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Langerin/CD207 is suitable for use in Immunofluorescence research applications.
Applications	IF
Host/Source	Mouse
Reactivity	Mouse/Rat/Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	310F7.02/HD26
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-546
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	