

Anti-Langerin/CD207 antibody [929F3.01] {Biotin} (STJ16102249)

STJ16102249

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Langerin/CD207 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rat
Reactivity	Human/Mouse/Rat/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	929F3.01
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence