

Anti-Langerin/CD207 antibody [306G9.01/HD24] {Alexa-fluor-647} (STJ16102258)

STJ16102258

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Langerin/CD207 is suitable for use in FACS research applications.
Applications	FACS
Host/Source	Mouse
Reactivity	Mouse/Rat/Human/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	306G9.01/HD24
Concentration	0.5 mg/mL
Conjugation	Alexa-fluor-647
Purification	Affinity Purified
Dilution Range	
Formulation	TRIS, PBS, pH7 and antibiotics
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence