

Anti-Langerin/CD207 antibody [306G9.01/HD24] {Biotin} (STJ16102259)
STJ16102259**GENERAL INFORMATION**

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Langerin/CD207 is suitable for use in FACS, Immunohistochemistry and Immunohistochemistry research applications. |
| Applications | FACS/IHC-F/IHC-P |
| Host/Source | Mouse |
| Reactivity | Mouse/Rat/Human/Pig |

PRODUCT PROPERTIES

| | |
|-----------------------|----------------------------------|
| Clonality | Monoclonal |
| Clone ID | 306G9.01/HD24 |
| Concentration | 0.5 mg/mL |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | |
| Formulation | TRIS, PBS, pH7 and antibiotics |
| Isotype | IgG1 |
| Storage | Store at -20°C for up to 1 year. |
| Instruction | |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence